

# LinguaGraph 3D Electropalatography System



3D Electropalatography (3D EPG) uses state-of-the-art computer animation and 3D modelling techniques to illustrate where the tongue touches the hard-palate during speech. The 3D model, which can be rotated 360°, provides direct articulatory information in real-time.

## How 3D Electropalatography works

3D Electropalatography determines tongue-palate contact by using a special artificial palate that contains 62 electrodes embedded on its lingual surface. The palate is custom-made for each user and simply clips to the upper teeth. The wires from the palate are connected to the EPG unit, which is worn around the neck.




A small electrical signal is passed to the user, which is conducted through the tongue to any touched electrodes. The detected tongue-palate contact is transferred to a computer and used to animate a 3D model of the tongue in real-time. The figure below illustrates six tongue-palate contact shapes observed in normal English-speakers when producing lingual consonants:




## 3D Electropalatography palate

3D Electropalatography is compatible with the Reading electropalatography palate. Each palate, sold separately, is custom-made to ensure a comfortable fit. Please contact us for details about our EPG palate manufacturing service. We also supply EPG palate kits that enable orthodontic technologists to manufacture palates locally.



## 3D software



LinguaView 3D incorporates state-of-the-art computer animation and 3D modelling techniques:

- Preview 3D tongue-palate contact in real-time
- Record 3D tongue-palate contact at up to 200 FPS with fully synchronised audio
- 3D tongue-palate contact template library
- 3D tongue-palate contact template editor
- 360° photorealistic and stylised views
- Upper and lower jaw with adjustable transparency
- File format compatible with icSpeech Professional Edition
- Compatible with Windows 64-bit 10 and 8



[www.rose-medical.com](http://www.rose-medical.com)



[enquiries@rose-medical.com](mailto:enquiries@rose-medical.com)



+44 (0) 1227 786 816

**RoseMedical**  
Speech Therapy Software & Instrumentation