

EAR, NOSE IRRIGATOR (Unit Attachable)

USER MANUAL (Ver.1)

USER MANUAL

□ MIR-2000



Ear, Nose Irrigator (Unit attachable)

1 Composition



- Saline Solution Tank
Volume Capacity : 1.2L ~ 5L



- Spray Gun



Ear tip



Nose tip



- Temperature Control
Setting and control temperature of saline solution.



- Spray Pressure Control
Control spray pressure

2 Description

* Power Consumption	220V, 50/60Hz, 300W
* Suitable Saline Temp	35°C (Max : 40°C)
* Saline Tank Volume	1.2L ~ 5L
* Spray Pressure	Max 1.5kgf/Cm
* Dimension	300(W) x 340(D) x 325(H)mm

3 Properties

This device made by basic medical knowledge that 35°C saline solution is optimum to washing of human tissue.

It consisted of saline tank which always keep 35°C temperature with safety tools.

Some people think that its running cost wold be high, due to saline price, but it doesn't matter because everyone can make saline solution with salt 5 ~ 10%. It is useful to clean up ear and nose part and/or other E.N.T. purposes.

4 How to attach to the Unit

a) Attach the arm set to the pole



b) Fix the irrigator shelf on the arm



c) Put the irrigator on the shelf safely



	 <p data-bbox="1053 1400 1420 1444">KAZAMAENT CO., LTD.</p>
---	--

<table border="1"><tr><td data-bbox="327 1635 406 1713">EC</td><td data-bbox="406 1635 486 1713">REP</td></tr></table>	EC	REP	
EC	REP		