

ENT TREATMENT FULL HD VISUAL SYSTEM : NET-260



1. Feature

1. Upgraded to high resolution of Full HD system.
2. Designed for easy and convenient diagnosis.
3. Very convenient gun type of camera for ENT diagnosis.
4. Freezing system enables image splitting and storing picture and video by USB port.
5. LAN port can be useful to be connected through distributor.
6. Easier lamp on/off and freezing by foot switch.

2. Specification.

1. Light Source (LED 60W)
2. Full HD Camera System
3. CCD 1/3" Camera (Auto Focus Control System)



Gun type



Straight type

3. Composition

1. Light Source 1
2. Color Monitor (Option) 1
3. Light Guide 1
4. CCD Camera (Kazama type) 1
5. Freezing System 1
6. USB Port (Max. 32G) 1
7. LAN Port 1
8. Foot Switch 1
9. Camera Adaptor 1
10. Operating & Service Manual 1



Front



Back

4. Picture Image

