



DIAGNOSIS VNG Ulmer

Evaluation by videonystagmography

The VNG Ulmer system enables to analyze the nystagmus at low and medium frequencies. It offers a large range of tests including vestibular and ocular-motor tests.

KEY BENEFITS

- Upgradable and customizable system
- Complete software with vestibular and ocular-motor tests
- Simple and intuitive
- Precise analyze of vestibular functions
- System designed by a doctor for doctors

The VNG Ulmer system is designed as a **modular system**. Each doctor can choose which tests to use so that the tool is **best suited for the practice** as well as to the associated camera (wired or wireless).

MODULES AVAILABLE:

- Nystagmus
- Caloric
- Saccade
- Smooth pursuit
- Gaze
- Optokinetic
- Kinetic

LATEST SOFTWARE VERSION FEATURES:

- **New interface** : intuitive and updated
- **New printout reports**
- **Easy customization of protocols** directly through software
- « **Plug and Play** » update : automatic back up of all old configurations during set up
- **Audio and video recording** with SmartCam camera (from 5.17.25 version)
- Available on **laptop** (new compatible cameras)

CONNECTIVITY:

- Available on Windows PC



DISPOSABLES:

- Goggles Flex with VNS4X wired camera
- Goggles Flex or Xpress Mask with SmartCam wireless camera

OPTIONS:

- **Binocular analysis** (with VISIO mask)
- **Rotational tests** (with MED4 chair)
- **Test of the otolith function** (with SVV module)
- **Detection of vibration-induced nystagmus** (with WIB3F vibrator)
- **Caloric tests** (with air or water irrigators)